

## BWB Product Comparison Information

[www.bwbtech.com](http://www.bwbtech.com)

	XP	XP Plus	BIO	BIO AV	Nuclear	Sugar	Lithium	Syn Fuels	Soil
<b>Sodium analysis</b>	✓	✓	✓	✓	✗	✓	✗	✓	✓
<b>Potassium analysis</b>	✓	✓	✓	✓	✗	✓	✗	✓	✓
<b>Lithium analysis</b>	✓	✓	✓	✓	✓	✗	✓	✓	✗
<b>Calcium analysis</b>	✓	✓	✓	✗	✗	✗	✗	✓	✓
<b>Barium analysis</b>	✓	✗	✗	✗	✗	✗	✗	✗	✗
<b>Lithium internal standard</b>	✗	✓	✓	✓	✗	✓	✗	✓	✓
<b>Cesium internal standard</b>	✗	✓	✓	✓	✓	✗	✓	✗	✗
<b>Simultaneous detection</b>	Na, K, Li, Ca, Ba	Na, K, Li, Ca	Na & K or Li & Ca	Na, K, Li	Li (across 4 ranges)	Na & K	Li (across 4 ranges)	Na, K, Li, Ca	Na, K, Ca
<b>Calibration points per Ion (excl blank)</b>	Up to 10	Up to 10	Up to 5	Up to 5	Up to 10 per channel	Up to 10	Up to 10 per channel	Up to 10	Up to 10
<b>Calibration curves (Saved methods)</b>	1	1	4	2	1	1	1	1	1
<b>Multi element calibration curve</b>	✓	✓	✓	✓	✗	✓	✗	✓	✓
<b>Internal results storage</b>	✓	✓	✓	✓	✓	✓	✓	✓	✓
<b>AFHS Link</b>	option	✓	✓	✓	option	option	option	option	option
<b>Integrated thermal printer</b>	option	✓	✓	✓	option	option	option	option	✓
<b>0 - 2.5v Chart output</b>	✓	✓	✓	✓	option	option	✓	✓	✓
<b>4 - 20mA Output</b>	option	option	option	option	✓	✓	option	option	option
<b>Collection cup</b>	option	option	option	option	✓	✓	✓	option	option
<b>Internal compressor</b>	✓	✓	✓	✓	✓	✓	✓	✓	✓
<b>Ready to use out of the box</b>	✓	✓	✓	✓	✓	✓	✓	✓	✓
<b>RS232 Datafeed output</b>	✗	✗	✗	✗	✗	✓	✗	✗	✗
<b>PC Software included</b>	✓	✓	✓	✓	✓	✓	✓	✓	✓
<b>PDF Report generation</b>	✓	✓	✓	✓	✓	✓	✓	✓	✓
<b>CSV report generation</b>	✓	✓	✓	✓	✓	✓	✓	✓	✓
<b>IQ OQ PQ possible</b>	✓	✓	✓	✓	✓	at request	✓	at request	at request
<b>IQ OQ PQ kit included</b>	✗	✗	✓	✓	✗	✗	✗	✗	✗